



ATC-04 CART STAND

IMPLANT/SURGERY SYSTEM CONSOLE



*** NOTE:**
BEFORE MOUNTING ASEPTICO
IMPLANT /SURGERY UNITS
AEU-17B AND AEU-20,
REMOVE RUBBER FEET THAT ARE
INSTALLED ON THE BOTTOM OF
THE CHASSIS.



8333 216th Street SE • Woodinville, WA 98072
1-800-426-5913 • Fax: 425-487-2608
www.aseptico.com

P/N: 420540 • Rev. A • ECO 11089 • 10/2005